


C-Fiber

Femtosecond Fiber Laser 1560 nm




Menlo Systems' fiber-based femtosecond laser sources integrate the latest achievements in fiber technology into easy-to-use products. Menlo Systems' unique figure 9[®] design results in reproducible and long-term stable operation. It is based on the well-established nonlinear optical loop mirror (NOLM) mode locking mechanism. Both oscillator and amplifier use polarization maintaining (PM) fiber components only, ensuring excellent stability and low-noise operation. The laser is maintenance free, user installed and ready to use at the press of a single button. Customize your laser with the available options to match the requirements of your application.

PERFORMANCE DATA


Amplitude noise

< 0.5% rms (over 24h)



Reproducibility

Identical and consistent laser performance



MenloSystems



KEY SPECIFICATIONS

- Wavelength 1560 nm
- Output Power >500 mW
- Pulse Length <90 fs
- Repetition Rate 50-250 MHz

APPLICATIONS

- Synchronization and Timing
- Ultrafast Spectroscopy
- Supercontinuum Generation
- Material Characterization
- Testing at Telecom Wavelengths

FEATURES

- High Stability
- Low Amplitude and Phase Noise
- All-PM Solution
- Single Mode-Lock State
- Menlo figure 9[®] Technology

OPTIONS

- **SYNC100**
Repetition Rate Synchronization
Tunable cavity length by high-bandwidth piezo-controlled synchronization
- **RRE-SYNCR0**
Repetition Rate Stabilization
Feedback electronics to phase lock pulses to an external clock (see separate data sheet for more details)
- **VARIO**
User-Defined Repetition Rate
Factory-set value selectable in the 50-250 MHz range
- **MULTIBRANCH**
Additional Seed Ports
Seeding of multiple amplifiers with optional subsequent frequency conversion to cover multiple wavelengths

C-Fiber

MenloSystems

Femtosecond Fiber Laser 1560 nm

| SPECIFICATIONS | C-FIBER | C-FIBER HIGH POWER |
|------------------------------------|--|---------------------------|
| Center Wavelength | 1560 nm ± 20 nm | 1560 nm ± 20 nm |
| Average Power | >100 mW | >500 mW |
| Pulse Energy | >1 nJ | >5 nJ |
| Pulse Width | <90 fs | |
| Repetition Rate | 100 MHz (50-250 MHz with VARIO)* | |
| Repetition Rate Instability | <1 ppm over 20 hours at constant temperature | |
| Timing Jitter | <2 fs [rms, 10 kHz.. 10 MHz] | |
| Output Port | fiber-coupled (FC/APC) | free space |
| Additional Fiber-Coupled Seed Port | 1 (up to 4 with MULTIBRANCH) | |
| Polarization | linear, PM fiber | linear, s-polarized |
| Beam Height | n.a. | 102 mm |

*Please inquire for your specific combinations of average power, pulse duration and repetition rate.

REQUIREMENTS

| | | |
|--------------------------------|---|---|
| Operating Voltage | 100/115/230 VAC | |
| Frequency | 50 to 60 Hz | |
| Power Consumption | 120 VA | |
| Cooling Requirements | no water cooling is required | |
| Laser Head Stabilization | actively temperature stabilized | |
| Operating Temperature | 15 °C - 35 °C | |
| Laser Head Dimensions/Weight | 413 x 90 x 178 mm ³ / 7 kg | 413 x 120 x 178 mm ³ / 10 kg |
| Control Unit Dimensions/Weight | 448 x 132 x 437 mm ³ / 10 kg | 448 x 132 x 437 mm ³ / 12 kg |
| Warm-Up Time | <60 s | |

ORDERING INFORMATION

| | | |
|---------------------|---------|--------------------|
| Product Code | C-Fiber | C-Fiber HIGH POWER |
|---------------------|---------|--------------------|

Please call for pricing. Specifications are subject to change without notice. Custom modifications are available, please inquire.

MenloSystems



Invisible laser radiation
avoid exposure to beam
Class 3B laser

Menlo Systems GmbH
T+49 89 189 166 0
sales@menlosystems.com

Menlo Systems, Inc.
T+1 973 300 4490
ussales@menlosystems.com

Thorlabs, Inc.
T+1 973 579 7227
sales@thorlabs.com



www.menlosystems.com

www.frequencycomb.com

D-CFiber-EN 05/05/17