

TERA Image

THz Imaging Extension



SPECIFICATIONS

Travel range	150 mm x 150 mm (ask us for longer ranges)	
Resolution x-direction	<500 μm^* (stage resolution <100 μm)	
Resolution y-direction	<500 μm^* (stage resolution <100 μm)	
Travel speed	TERA K15	up to 20 data pixels/second
	TeraSmart	
	TERA ASOPS	up to 100 data pixels/second (ask for fast imaging option!)
	TERA OSCAT	

*image resolution is limited by the size of the THz focal diameter

DIMENSIONS AND WEIGHT

Width x Height	380 mm x 430 mm
Weight	8 kg

SYSTEM COMPONENTS

	XY-scanning unit with controllers
	LabView-based measurement software
	EPINA ImageLab software for image reconstruction and analysis

ORDERING INFORMATION

Product Code	TERA Image
--------------	------------



Menlo Systems GmbH
T+49 89 189 166 0
sales@menlosystems.com

Menlo Systems, Inc.
T+1 973 300 4490
ussales@menlosystems.com

Thorlabs, Inc.
T+1 973 579 7227
sales@thorlabs.com

